## **Technical Data Sheet**

# Aquaperm® AP 9680 Permanent Marker Ink



## Water Based White Permanent Marker Ink

#### **General Information**

Aquaperm® AP 9680 is the world's first white ink that combines the high covering strength of a valve action marker ink and the easy handling of a reservoir marker. It has a good adhesion on many surfaces and excellent light-fastness properties.

#### **Suitable Surfaces**

For all kind of smooth and porous surfaces.

**Applicable Standards & Regulations on Request** 

## **Physical Data**

Туре	Colour	Viscosity at 20°C mPa s	Density at 20°C g/cm3 ± 0.2	pH-Value
Aquaperm® AP 9680	white	30.0 – 60.0	1.50	7.5 – 9.0

#### **Technical Advice**

**Suitable Components** 

Based on our tests we recommend following components:

Nibs: AuBex MO454 Reservoirs: Essentra

a) XPE 438156 (fibre density: 0.10 g/cm³) for maximum ink lay down, and acceptable sedimentation.

b) XPE 438349 (fibre density:0.15 g/cm³) for minimal sedimentation and an acceptable ink lay down.

Barrels, Caps, Plugs: polypropylene

If used in valve action markers, a corrosion safe ball (SS or glass) inside the barrel grants a good distribution of the pigments by shaking the marker before use.

It is mandatory to approve the compatibility of the ink and the components as well as the performance of the writing system.

### Storage, Handling & Transportation

Before use of our product the filling station has to be maintained and cleaned carefully. The inks are sensitive to any kind of contamination. Cross-contamination of the inks may lead to a colour mismatch or will even lead to a defective marker or shorter shelf life.

Store our product frost protected and avoid direct sunlight. The recommended storage temperature is between 10°C and 30°C. Storage temperatures of above 30°C for several days bear the risk of a reduced shelf life of the ink.

The drums should be closed tightly during storage. The shelf life of the ink in original sealed containers is one year.

Stir well before use and stir ink during ink injection in order to avoid pigment sedimentation during filling process. Markers should be stored horizontally.

#### **Packaging**

standard packaging: 20kg plastic canister (further packaging units on request).

## **Colour Chart**



## Produced by

DOKUMENTAL GmbH & Co KG Wöllnerstraße 26 D-67065 Ludwigshafen

Phone: +49(0)621/37702-200
Fax: +49(0)621/37702-391
Mail: info@dokumental.de
www.dokumental.de